

Product	Description	EUR
NexiusZoom - Binocular models		
NZ.1902-P	Binocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with pillar stand. Incident and transmitted 3 W LED illuminations	747,00 €
	Pair wide field eyepieces HWF10x/22mm without eyecups Binocular head with 45° inclined tubes Diopter adjustment ± 5° on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Ergonomically pillar stand. Maximum object height 148 mm Supplied with opaque stage plate, white/black stage plate and 2 object clamps Incident and transmitted 3W LED illumination with internal power supply Both illumination intensities can be adjusted separately Paint with Antibacterial Protection Layer Supplied with power cord, dust cover and user manual in Styrofoam case	
NZ.1902-PG	Binocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with pillar stand. Two incident LED illuminations with Gooseneck and one transmitted 3 W LED illumination	836,00 €
	Pair wide field eyepieces HWF10x/22mm without eyecups Binocular head with 45° inclined tubes Diopter adjustment ± 5° on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Ergonomically pillar stand. Maximum object height 148 mm Supplied with opaque stage plate, white/black stage plate and 2 object clamps Two incident LED illuminations with Gooseneck and one transmitted 3 W LED illumination. Internal power supply Both illumination intensities can be adjusted separately Paint with Antibacterial Protection Layer Supplied with power cord, dust cover and user manual in Styrofoam case	
NZ.1902-S	Binocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with rack and pinion stand. Incident and transmitted 3 W LED illuminations	760,00 €
	Pair wide field eyepieces HWF10x/22mm without eyecups Binocular head with 45° inclined tubes Diopter adjustment ± 5° on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Ergonomically rack and pinion stand. Maximum object height 58 mm Supplied with opaque stage plate, white/black stage plate and 2 object clamps Incident and transmitted 3W LED illumination with internal power supply Both illumination intensities can be adjusted separately Paint with Antibacterial Protection Layer Supplied with power cord, dust cover and user manual in Styrofoam case	

Product	Description	EUR
NZ.1902-U	Binocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with universal one-arm stand. Without illumination	868,00 €
	Pair wide field eyepieces HWF10x/22mm without eyecups Binocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Universal single-arm stand. Maximum object height 213 mm Without illumination Paint with Antibacterial Protection Layer Supplied with dust cover and manual	
NZ.1902-B	Binocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with double-arm boom stand. Without illumination	915,00 €
	Pair wide field eyepieces HWF10x/22mm without eyecups Binocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Double-arm boom stand. Maximum object height 213 mm Without illumination Paint with Antibacterial Protection Layer	
NZ.1902-A	Binocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with articulated stand for table mounting. Without illumination	1.038,00 €
	Pair wide field eyepieces HWF10x/22mm without eyecups Binocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Articulated stand for table mounting without illumination Paint with Antibacterial Protection Layer	WSL

Product	Description	EUR
NZ.1902-M	Binocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with pillar stand. Incident 3 W LED illumination and transmitted 3 W LED illuminations with rotating mirror	1.054,00 €
	Pair wide field eyepieces HWF10x/22mm without eyecups Binocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Ergonomically pillar stand with rotating mirror. Maximum object height 148 mm Supplied with opaque stage plate, white/black stage plate and 2 object clamps Incident 3W LED illumination transmitted 3W LED illumination with rotating mirror. Both with internal power supply Both illumination intensities can be adjusted separately Supplied with power cord, dust cover and user manual in Styrofoam case	
NexiusZoom - Trinocular models		
NZ.1903-P	Trinocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with pillar. Incident and transmitted 3 W LED illuminations	805,00 €
	Pair wide field eyepieces HWF10x/22mm without eyecups Trinocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Ergonomically pillar stand. Maximum object height 148 mm Supplied with opaque stage plate, white/black stage plate and 2 object clamps Incident and transmitted 3W LED illumination with internal power supply. Both illumination intensities can be adjusted separately Paint with Antibacterial Protection Layer Supplied with dust cover in Styrofoam case	
NZ.1903-PG	Trinocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with pillar stand. Two incident LED illuminations with Gooseneck and one transmitted 3 W LED illumination	894,00 €
	Pair wide field eyepieces HWF10x/22mm without eyecups Trinocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Ergonomically pillar stand. Maximum object height 148 mm Supplied with opaque stage plate, white/black stage plate and 2 object clamps Two incident LED illuminations with Gooseneck and one transmitted 3 W LED illumination. Internal power supply Both illumination intensities can be adjusted separately Paint with Antibacterial Protection Layer Supplied with power cord, dust cover and user manual in Styrofoam case	

Product	Description	EUR
NZ.1903-S	Trinocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with rack and pinion stand. Incident and transmitted 3 W LED illuminations	817,00 €
	Pair wide field eyepieces HWF10x/22mm without eyecups Trinocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Ergonomically rack and pinion stand. Maximum object height 58 mm Supplied with opaque stage plate, white/black stage plate and 2 object clamps Incident and transmitted 3W LED illumination with internal power supply Both illumination intensities can be adjusted separately Paint with Antibacterial Protection Layer Supplied with dust cover in Styrofoam case	
NZ.1903-U	Trinocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with universal one-arm stand. Without illumination	924,00 €
	Pair wide field eyepieces HWF10x/22mm without eyecups Trinocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Universal single-arm stand. Maximum object height 213 mm Without illumination Paint with Antibacterial Protection Layer	
NZ.1903-B	Trinocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with double-arm boom stand. Without illumination	973,00 €
	Pair wide field eyepieces HWF10x/22mm without eyecups Trinocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Double-arm boom stand. Maximum object height 213 mm Without illumination Paint with Antibacterial Protection Layer	
NZ.1903-A	Trinocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with articulated stand for table mounting. Without illumination	1.098,00 € WSL
	Pair wide field eyepieces HWF10x/22mm without eyecups Trinocular head with 45° inclined tubes Diopter adjustment on $\pm 5^\circ$ two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Articulated stand for table mounting. Without illumination Paint with Antibacterial Protection Layer	

Product	Description	EUR
NZ.1903-AP	Trinocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with articulated stand with heavy tabletop stand. Without illumination	1.121,00 €
	Pair wide field eyepieces HWF10x/22mm without eyecups Trinocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Articulated stand with heavy tabletop stand. Without illumination	
NZ.1903-M	Trinocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with pillar stand. Incident 3 W LED illumination and transmitted 3 W LED illuminations with rotating mirror	1.111,00 €
	Pair wide field eyepieces HWF10x/22mm without eyecups Trinocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Ergonomically pillar stand with rotating mirror. Maximum object height 148 mm Supplied with opaque stage plate, white/black stage plate and 2 object clamps Incident 3W LED illumination transmitted 3W LED illumination with rotating mirror. Both with internal power supply Both illumination intensities can be adjusted separately Paint with Antibacterial Protection Layer Supplied with power cord, dust cover and user manual in Styrofoam case	
NexiusZoom ESD - Binocular models		
NZ.1902-P-ESD	Binocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with pillar. Incident and transmitted 3 W LED illuminations. Antistatic (ESD) version	958,00 €
	Pair wide field eyepieces HWF10x/22mm without eyecups Binocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Ergonomically pillar stand. Maximum object height 148 mm Supplied with opaque stage plate, white/black stage plate and 2 object clamps Incident and transmitted 3W LED illumination with internal power supply Both illumination intensities can be adjusted separately The body and stand are coated with a static dissipative paint eliminating harmful electrostatic discharges Supplied with dust cover in Styrofoam case	

Product	Description	EUR
NZ.1902-S-ESD	Binocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with rack and pinion stand. Incident and transmitted 3 W LED illuminations. Antistatic (ESD) version	1.007,00 €
	Pair wide field eyepieces HWF10x/22mm without eyecups Binocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Ergonomically rack and pinion stand. Maximum object height 58 mm Supplied with opaque stage plate, white/black stage plate and 2 object clamps Incident and transmitted 3W LED illumination with internal power supply Both illumination intensities can be adjusted separately The body and stand are coated with a static dissipative paint eliminating harmful electrostatic discharges Supplied with dust cover in Styrofoam case	
NZ.1902-U-ESD	Binocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with universal one-arm stand. Without illumination. Antistatic (ESD) version	1.226,00 €
	Pair wide field eyepieces HWF10x/22mm without eyecups Binocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Universal single-arm stand. Maximum object height 213 mm Without illumination The body and stand are coated with a static dissipative paint eliminating harmful electrostatic discharges Supplied with dust cover and manual	
NZ.1902-B-ESD	Binocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with double-arm boom stand. Without illumination. Antistatic (ESD) version	1.364,00 €
	Pair wide field eyepieces HWF10x/22mm without eyecups Binocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Double-arm boom stand. Maximum object height 213 mm Without illumination The body and stand are coated with a static dissipative paint eliminating harmful electrostatic discharges Supplied with dust cover and manual	

Product	Description	EUR
NexiusZoom ESD - Trinocular models		
NZ.1903-P-ESD	Trinocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with pillar. Incident and transmitted 3 W LED illuminations. Antistatic (ESD) version	1.040,00 €
	Pair wide field eyepieces HWF10x/22mm without eyecups Trinocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Ergonomically pillar stand. Maximum object height 148 mm Supplied with opaque stage plate, white/black stage plate and 2 object clamps The body and stand are coated with a static dissipative paint eliminating harmful electrostatic discharges Incident and transmitted 3W LED illumination with internal power supply Both illumination intensities can be adjusted separately Supplied with dust cover in Styrofoam case	
NZ.1903-S-ESD	Trinocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with rack and pinion stand. Incident and transmitted 3 W LED illuminations. Antistatic (ESD) version	1.090,00 €
	Pair wide field eyepieces HWF10x/22mm without eyecups Trinocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Ergonomically rack and pinion stand. Maximum object height 58 mm Supplied with opaque stage plate, white/black stage plate and 2 object clamps The body and stand are coated with a static dissipative paint eliminating harmful electrostatic discharges Incident and transmitted 3W LED illumination with internal power supply Both illumination intensities can be adjusted separately Supplied with dust cover in Styrofoam case	
NZ.1903-U-ESD	Trinocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with universal one-arm stand. Without illumination. Antistatic (ESD) version	1.310,00 €
	Pair wide field eyepieces HWF10x/22mm without eyecups Trinocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Universal one-arm stand. Maximum object height 213 mm Without illumination The body and stand are coated with a static dissipative paint eliminating harmful electrostatic discharges	

Product	Description	EUR
NZ.1903-B-ESD	Trinocular stereo zoom microscope NexusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with double-arm boom stand. Without illumination. Antistatic (ESD) version	1.447,00 €
	Pair wide field eyepieces HWF10x/22mm without eyecups Trinocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Double-arm boom stand. Maximum object height 213 mm Without illumination The body and stand are coated with a static dissipative paint eliminating harmful electrostatic discharges	
	NexusZoom EVO - Binocular models	
NZ.1702-P	Binocular stereo zoom microscope NexusZoom EVO, 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with pillar. Incident and transmitted 3 W LED illuminations	950,00 €
	Pair wide field eyepieces HWF10x/23mm without eyecups Binocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with click stops Field of view from 35 mm to 4.2 mm. Working distance 110 mm Ergonomically pillar stand. Maximum object height 148 mm Supplied with opaque stage plate, white/black stage plate and 2 object clamps Incident and transmitted 3W LED illumination with internal power supply Both illumination intensities can be adjusted separately Paint with Antibacterial Protection Layer Supplied with power cord, dust cover and user manual in Styrofoam case	
NZ.1702-PG	Binocular stereo zoom microscope NexusZoom EVO, 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with pillar. Two incident LED illuminations with Gooseneck and one transmitted 3 W LED illumination	1.052,00 €
	Pair wide field eyepieces HWF10x/23mm without eyecups Binocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with click stops Field of view from 35 mm to 4.2 mm. Working distance 110 mm Ergonomically pillar stand. Maximum object height 148 mm Supplied with opaque stage plate, white/black stage plate and 2 object clamps Two incident LED illuminations with Gooseneck and one transmitted 3 W LED illumination. Internal power supply Both illumination intensities can be adjusted separately Paint with Antibacterial Protection Layer Supplied with power cord, dust cover and user manual in Styrofoam case	

Product	Description	EUR
NZ.1702-S	Binocular stereo zoom microscope NexiusZoom EVO, 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with rack and pinion stand. Incident and transmitted 3 W LED illuminations	981,00 €
	Pair wide field eyepieces HWF10x/23mm without eyecups Binocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with click stops Field of view from 35 mm to 4.2 mm. Working distance 110 mm Ergonomically rack and pinion stand. Maximum object height 58 mm Supplied with opaque stage plate, white/black stage plate and 2 object clamps Incident and transmitted 3W LED illumination with internal power supply Both illumination intensities can be adjusted separately Paint with Antibacterial Protection Layer Supplied with power cord, dust cover and user manual in Styrofoam case	
NZ.1702-U	Binocular stereo zoom microscope NexiusZoom EVO, 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with universal one-arm stand. Without illumination	1.087,00 €
	Pair wide field eyepieces HWF10x/23mm without eyecups Binocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with click stop Field of view from 33 mm to 4.9 mm. Working distance 110 mm Universal single-arm stand. Maximum object height 213 mm Without illumination Paint with Antibacterial Protection Layer Supplied with dust cover and manual	
NZ.1702-B	Binocular stereo zoom microscope NexiusZoom EVO, 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with double-arm boom stand. Without illumination	1.144,00 €
	Pair wide field eyepieces HWF10x/23mm without eyecups Binocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with click stops Field of view from 35 mm to 4.2 mm. Working distance 110 mm Double-arm boom stand. Maximum object height 213 mm Without illumination Paint with Antibacterial Protection Layer	

Product	Description	EUR	
NZ.1702-A	Binocular stereo zoom microscope NexiusZoom EVO, 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with articulated stand for table mounting. Without illumination	1.211,00 €	WSL
	Pair wide field eyepieces HWF10x/23mm without eyecups		
	Binocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with click stops Field of view from 35 mm to 4.2 mm. Working distance 110 mm Articulated stand for table mounting. Without illumination		
NZ.1702-AP	Binocular stereo zoom microscope NexiusZoom EVO, 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with articulated stand with heavy stand. Without illumination	1.234,00 €	WSL
	Pair wide field eyepieces HWF10x/23mm without eyecups		
	Binocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with click stops Field of view from 35 mm to 4.2 mm. Working distance 110 mm Articulated stand with heavy stand. Without illumination Paint with Antibacterial Protection Layer		
NZ.1702-M	Binocular stereo zoom microscope NexiusZoom EVO, 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with pillar. Incident 3 W LED illumination and transmitted 3 W LED illuminations with rotating mirror	1.519,00 €	
	Pair wide field eyepieces HWF10x/23mm without eyecups		
	Binocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with click stops Field of view from 35 mm to 4.2 mm. Working distance 110 mm Ergonomically pillar stand. Maximum object height 148 mm Supplied with opaque stage plate, white/black stage plate and 2 object clamps Incident 3W LED illumination transmitted 3W LED illumination with rotating mirror. Both with internal power supply Both illumination intensities can be adjusted separately Paint with Antibacterial Protection Layer Supplied with power cord, dust cover and user manual in Styrofoam case		

Product	Description	EUR
NexiusZoom EVO - Trinocular models		
NZ.1703-P	Trinocular stereo zoom microscope NexiusZoom EVO, 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with pillar. Incident and transmitted 3 W LED illuminations	1.007,00 €
	Pair wide field eyepieces HWF10x/23mm without eyecups Trinocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with click-stops Field of view from 35 mm to 4.2 mm. Working distance 110 mm Ergonomically pillar stand. Maximum object height 148 mm Supplied with opaque stage plate, white/black stage plate and 2 object clamps Incident and transmitted 3W LED illumination with internal power supply. Both illumination intensities can be adjusted separately Paint with Antibacterial Protection Layer Supplied with dust cover in Styrofoam case	
NZ.1703-PG	Trinocular stereo zoom microscope NexiusZoom EVO, 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with pillar. Two incident LED illuminations with Gooseneck and one transmitted 3 W LED illumination	1.180,00 €
	Pair wide field eyepieces HWF10x/23mm without eyecups Trinocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with click-stops Field of view from 35 mm to 4.2 mm. Working distance 110 mm Ergonomically pillar stand. Maximum object height 148 mm Supplied with opaque stage plate, white/black stage plate and 2 object clamps Two incident LED illuminations with Gooseneck and one transmitted 3 W LED illumination. Internal power supply Both illumination intensities can be adjusted separately Paint with Antibacterial Protection Layer Supplied with dust cover in Styrofoam case	
NZ.1703-S	Trinocular stereo zoom microscope NexiusZoom EVO, 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with rack and pinion stand. Incident and transmitted 3 W LED illuminations	1.038,00 €
	Pair wide field eyepieces HWF10x/23mm without eyecups Trinocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with click-stops Field of view from 35 mm to 4.2 mm. Working distance 110 mm Ergonomically rack and pinion stand. Maximum object height 58 mm Supplied with opaque stage plate, white/black stage plate and 2 object clamps Incident and transmitted 3W LED illumination with internal power supply Both illumination intensities can be adjusted separately Paint with Antibacterial Protection Layer Supplied with dust cover in Styrofoam case	

Product	Description	EUR
NZ.1703-U	Trinocular stereo zoom microscope NexiusZoom EVO, 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with universal one-arm stand. Without illumination	1.144,00 €
	Pair wide field eyepieces HWF10x/23mm without eyecups Trinocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with click-stops Field of view from 35 mm to 4.2 mm. Working distance 110 mm Universal single-arm stand. Maximum object height 213 mm Without illumination Paint with Antibacterial Protection Layer	
NZ.1703-B	Trinocular stereo zoom microscope NexiusZoom EVO, 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with double-arm boom stand. Without illumination	1.202,00 €
	Pair wide field eyepieces HWF10x/23mm without eyecups Trinocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with click-stops Field of view from 35 mm to 4.2 mm. Working distance 110 mm Double-arm boom stand. Maximum object height 213 mm Without illumination Paint with Antibacterial Protection Layer	
NZ.1703-A	Trinocular stereo zoom microscope NexiusZoom EVO, 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with articulated stand for table mounting. Without illumination	1.263,00 € WSL
	Pair wide field eyepieces HWF10x/23mm without eyecups Trinocular head with 45° inclined tubes Diopter adjustment on $\pm 5^\circ$ two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with click-stops Field of view from 35 mm to 4.2 mm. Working distance 110 mm Articulated stand without illumination	
NZ.1703-AP	Trinocular stereo zoom microscope NexiusZoom EVO, 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with articulated stand with heavy stand. Without illumination	1.285,00 € WSL
	Pair wide field eyepieces HWF10x/23mm without eyecups Trinocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with click-stops Field of view from 35 mm to 4.2 mm. Working distance 110 mm Articulated stand with heavy stand. Without illumination Paint with Antibacterial Protection Layer	

Product	Description	EUR
NZ.1703-M	Trinocular stereo zoom microscope NexiusZoom EVO, 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with pillar. Incident 3 W LED illumination and transmitted 3 W LED illuminations with rotating mirror	1.575,00 €
	Pair wide field eyepieces HWF10x/23mm without eyecups Trinocular head with 45° inclined tubes Diopter adjustment ± 5° on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with click-stops Field of view from 35 mm to 4.2 mm. Working distance 110 mm Ergonomically pillar stand. Maximum object height 148 mm Supplied with opaque stage plate, white/black stage plate and 2 object clamps Incident 3W LED illumination transmitted 3W LED illumination with rotating mirror. Both with internal power supply Both illumination intensities can be adjusted separately Paint with Antibacterial Protection Layer Supplied with power cord, dust cover and user manual in Styrofoam case	
	NexiusZoom for Gemology	
NZ.1902-GEML	Binocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with stand for gemology with 0-45° tilting arm for dark- and brightfield contrast, Built-in 100-240 V power supply, a 30W 6V halogen transmitted illumination and 1W LED incident illumination	1.820,00 € BTO
	Pair wide field eyepieces HWF10x/22mm without eyecups Binocular head with 45° inclined tubes Diopter adjustment ± 5° on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Stand for gemology with 0-45° tilting arm for dark- and brightfield contrast. Built-in 100-240 V power supply 30W 6V halogen transmitted illumination 1W LED incident illumination	

Product	Description	EUR
NZ.1902-GEMF	Binocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with stand for gemology with 0-45° tilting arm for dark- and brightfield contrast, Built-in 100-240 V power supply, a 30W 6V halogen transmitted illumination and a 7W fluorescent gas-discharge lamp incident illumination	1.740,00 € BTO
	Pair wide field eyepieces HWF10x/22mm without eyecups Binocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Stand for gemology with 0-45° tilting arm for dark- and brightfield contrast. Built-in 100-240 V power supply 30W 6V halogen transmitted illumination 7W fluorescent incident illumination	
NZ.1903-GEML	Trinocular stereo zoom microscope NexiusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with stand for gemology with 0-45° tilting arm for dark- and brightfield contrast, Built-in 100-240 V power supply, a 30W 6V halogen transmitted illumination and 1W LED incident illumination	1.979,00 € BTO
	Pair wide field eyepieces HWF10x/22mm without eyecups Trinocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Stand for gemology with 0-45° tilting arm for dark- and brightfield contrast. Built-in 100-240 V power supply 30W 6V halogen transmitted illumination 1W LED incident illumination	

Product	Description	EUR
NZ.1903-GEMF	Trinocular stereo zoom microscope NexusZoom, 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x with stand for gemology with 0-45° tilting arm for dark- and brightfield contrast, Built-in 100-240 V power supply, a 30W 6V halogen transmitted illumination and a 7W fluorescent gas-discharge lamp incident illumination	1.918,00 € BTO
	Pair wide field eyepieces HWF10x/22mm without eyecups Trinocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.67x to 4.5x zoom objective, magnification from 6.7x to 45x Field of view from 33 mm to 4.9 mm. Working distance 110 mm Stand for gemology with 0-45° tilting arm for dark- and brightfield contrast. Built-in 100-240 V power supply 30W 6V halogen transmitted illumination 7W fluorescent incident illumination Paint with Antibacterial Protection Layer	
NZ.1702-GEML	Binocular stereo zoom microscope NexusZoom EVO, 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with stand for gemology with 0-45° tilting arm for dark- and brightfield contrast, Built-in 100-240 V power supply, a 30W 6V halogen transmitted illumination and 1W LED incident illumination	2.000,00 €
	Pair wide field eyepieces HWF10x/23mm without eyecups Binocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with click-stops Field of view from 35 mm to 4.2 mm. Working distance 110 mm Stand for gemology with 0-45° tilting arm for dark- and brightfield contrast. Built-in 100-240 V power supply 30W 6V halogen transmitted illumination 1W LED incident illumination	
NZ.1702-GEMF	Binocular stereo zoom microscope NexusZoom EVO, 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with stand for gemology with 0-45° tilting arm for dark- and brightfield contrast, Built-in 100-240 V power supply, a 30W 6V halogen transmitted illumination and a 7W fluorescent gas-discharge lamp incident illumination	1.980,00 €
	Pair wide field eyepieces HWF10x/23mm without eyecups Binocular head with 45° inclined tubes Diopter adjustment $\pm 5^\circ$ on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with click-stops Field of view from 35 mm to 4.2 mm. Working distance 110 mm Stand for gemology with 0-45° tilting arm for dark- and brightfield contrast. Built-in 100-240 V power supply 30W 6V halogen transmitted illumination 7W fluorescent incident illumination	

Product	Description	EUR
NZ.1703-GEML	Trinocular stereo zoom microscope NexiusZoom EVO, 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with stand for gemology with 0-45° tilting arm for dark- and brightfield contrast, Built-in 100-240 V power supply, a 30W 6V halogen transmitted illumination and 1W LED incident illumination	2.186,00 €
	Pair wide field eyepieces HWF10x/23mm without eyecups	
	Trinocular head with 45° inclined tubes Diopter adjustment ± 5° on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with click-stops Field of view from 35 mm to 4.2 mm. Working distance 110 mm Stand for gemology with 0-45° tilting arm for dark- and brightfield contrast. Built-in 100-240 V power supply 30W 6V halogen transmitted illumination 1W LED incident illumination	
NZ.1703-GEMF	Trinocular stereo zoom microscope NexiusZoom EVO, 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with stand for gemology with 0-45° tilting arm for dark- and brightfield contrast, Built-in 100-240 V power supply, a 30W 6V halogen transmitted illumination and a 7W fluorescent gas-discharge lamp incident illumination	2.168,00 €
	Pair wide field eyepieces HWF10x/23mm without eyecups	
	Trinocular head with 45° inclined tubes Diopter adjustment ± 5° on two tubes. Interpupillary distance adjustable between 54 and 75 mm 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with click-stops Field of view from 35 mm to 4.2 mm. Working distance 110 mm Stand for gemology with 0-45° tilting arm for dark- and brightfield contrast. Built-in 100-240 V power supply 30W 6V halogen transmitted illumination 7W fluorescent incident illumination Paint with Antibacterial Protection Layer	
Heads for NexiusZoom		
NZ.5302	NexiusZoom binocular head 0.67-4.5x without head holder	562,00 €
NZ.5303	NexiusZoom trinocular head 0.67-4.5x without head holder	630,00 €
NZ.5312	NexiusZoom EVO binocular head 0.65-5.5x without head holder	755,00 €
NZ.5313	NexiusZoom EVO trinocular head 0.65-5.5x without head holder	812,00 €
Eyepieces for NexiusZoom		
NZ.6010	Pair of HWF 10x/22 mm eyepieces for NexiusZoom	171,00 €
NZ.6010-C	HWF 10x/22 mm eyepiece only with crosshair	101,00 €
NZ.6010-CM	HWF 10x/22 mm eyepiece with 10/100 micrometer and cross hair	114,00 €
NZ.6015	Pair of HWF 15x/16 mm eyepieces for NexiusZoom	182,00 €
NZ.6020	Pair of HWF 20x/12 mm eyepieces for NexiusZoom	182,00 €
NZ.6020-CM	HWF 20x/12 mm eyepieces with 10/100 micrometer and cross hair for NexiusZoom	194,00 €
NZ.6099	Pair of eyecups for NexiusZoom	6,10 €
NZ.6210	Pair of HWF 10x/23 mm eyepieces for NexiusZoom EVO	179,00 €
NZ.6210-C	HWF 10x/23 mm eyepiece only with crosshair for NexiusZoom EVO	101,00 €
NZ.6210-CM	HWF 10x/23 mm eyepiece only with 10/100 micrometer and cross hair for NexiusZoom EVO	114,00 €

Product	Description	EUR
Additional lenses for NexiusZoom		
NZ.8903	Additional 0.3x lens for NexiusZoom. Working distance 287 mm +/- 10 mm. Only suitable for models for A, AP, U, B, BC and M stands	132,00 €
NZ.8904	Additional 0.4x lens for NexiusZoom. Working distance 220 mm. Only suitable for P, PG, A, AP, U, B, BC and M stands	132,00 €
NZ.8905	Additional 0.5x lens for NexiusZoom. Working distance 183 mm. Only suitable for P, PG, A, AP, U, B, BC and M stands	132,00 €
NZ.8907	Additional 0.75x lens for NexiusZoom. Working distance 105 mm. Only suitable for P, PG, A, AP, U, B, BC and M stands	132,00 €
NZ.8915	Additional 1.5x lens for NexiusZoom. Working distance 53 mm	132,00 €
NZ.8920	Additional 2.0 lens for NexiusZoom. Working distance 34 mm	153,00 €
NZ.8950	Protection glass window for NexiusZoom	24,00 €
Stands for NexiusZoom		
NZ.9000	Ergo pillar stand without illumination	165,00 €
		
NZ.9005	Ergo pillar stand with adjustable transmitted 3W LED illumination	333,00 €
		
NZ.9010	Ergo rack & pinion stand with incident and transmitted 3W LED illuminators	322,00 €
		
NZ.9015	Ergo pillar stand with incident and transmitted 3W LED illuminators and head support	260,00 €
		
NZ.9018	Ergonomically designed pillar stand with two 3W goose-neck type incident LED illuminations on each side and a 3W transmitted LED illumination	342,00 €
		

Product	Description	EUR	
NZ.9020	Universal one-arm stand without NexiusZoom head holder	327,00 €	
			
NZ.9025	Black universal stand with <u>table clamp</u> without NexiusZoom head holder	450,00 €	WSL
	Black universal stand with table clamp without NexiusZoom head holder. Can be equipped with a 76 mm head holder for NexiusZoom (NZ.9081) or DZ series (65.981), 83 mm head holder (65.982) or an 85 mm head holder for ZE series (65.983)		
NZ.9027	Black universal stand with <u>heavy baseplate</u> without NexiusZoom head holder	475,00 €	WSL
	Black universal stand with heavy baseplate without NexiusZoom head holder. Can be equipped with a 76 mm head holder for NexiusZoom (NZ.9081) or DZ series (65.981), 83 mm head holder (65.982) or an 85 mm head holder for ZE series (65.983)		
NZ.9030	Universal double-arm boom stand without NexiusZoom head holder	381,00 €	
			
NZ.9032	Universal two-arm stand with table clamp without NexiusZoom holder	343,00 €	
			
NZ.9081	Head 76 mm holder for NexiusZoom and black universal stand NZ.9025 / NZ.9027	141,00 €	
NZ.9040	Darkfield attachment for NexiusZoom	91,00 €	
NZ.9042	Stand with rotating mirror and indirect transmitted LED illumination	717,00 €	
			
NZ.9045	Fall protection collar for NexiusZoom pillars	17,90 €	
NZ.9047	Articulated stand with adjustable tilt, swivel and rotating functions. With desktop clamp for mounting on 10-85 mm tabletops. Can also be fixed on tables with drilled holes. Supplied with 76 mm head holder for NexiusZoom stereo heads	331,00 €	
			
NZ.9090	NexiusZoom head holder for NZ.9000/9005/9015/9017/9020/9030 pillar stands	97,00 €	
NZ.9095	NexiusZoom head holder with fine coarse for NZ.9000/9005/9015/9017/9020/9030 pillar stands	195,00 €	
NZ.9095-ADD	NexiusZoom head holder with fine coarse ADD-ON for NZ.9000 / 9005 / 9015 / 9017 / 9020 / 9030 pillar stands. Only available with new microscopes	66,00 €	

Product	Description	EUR
Miscellaneous		
NZ.4300	Aluminium case for NexiusZoom microscopes with -P/-S/-PG stands	107,00 €
NZ.9505	Mechanical 180 x 155 mm X-Y stage with 75 x 75 mm translation and transparent glass plate. For mounting hole with diameter 95 mm. Only available with new microscope	189,00 €
NZ.9520	Polarization kit for NexiusZoom: 360° rotatable round stage with built-in polarization filter (NZ.9524) + analyzer in mount to be screwed under head (NZ.9525). For mounting hole with diameter 95 mm. Only available with new microscope	208,00 €
NZ.9524	360° rotatable round stage with built-in polarization filter for NexiusZoom	131,00 €
NZ.9525	360° rotatable analyzer in mount to be screwed under head of NexiusZoom	80,00 €
NZ.9570	Pair object clamps for standard stage	6,30 €
NZ.9572	Adjustable GEM object clamp for P, PL and S stands of NexiusZoom	41,00 €
NZ.9850	C-mount adapter with 0.5x lens for 1/2 inch cameras	171,00 €
NZ.9950	Standard Plexiglas-glass object plate, opaque, with black border	15,20 €
NZ.9956	Black/white object plate	30,00 €
NZ.9958	Standard glass object plate, opaque	15,20 €
NZ.9983-R	3W LED replacement unit for NexiusZoom, incident illumination	26,00 €
NZ.9983	3 W LED replacement unit for NexiusZoom, transmitted illumination	26,00 €
AE.9916	Euromex dust cover extra-large. For BioBlue (trino), EcoBlue (trino), BioBlue.Lab, Oxion, bScope, iScope, StereoBlue and NexiusZoom	10,30 €
AE.5130	Universal SLR adapter with built-in 2x lens for standard 23.2 mm tube. Needs T2 adapter	261,00 €
AE.5025	T2 ring for Nikon D SLR digital camera	38,00 €
AE.5040	T2 ring for Canon EOS SLR digital camera	38,00 €
AE.5168-NZ	X-Y stage with heating elements and controller for NZ microscope stage. Only with the purchase of a microscope	1.235,00 €
AE.5168-UNI	Plate with heating elements and controller, to be placed on top. Dimensions: 150 (W) x 130 (D) x 15 H) mm. With hole 40 x 20 mm	495,00 €
AE.1112	Micrometer 76 x 26 mm slide, 50 mm/500 divisions, 100 µm per division	76,00 €
PB.5245	Lens cleaning paper, 100 sheets per pack	3,20 €
PB.5274	Iso propyl alcohol 99%, 200 ml	29,00 €
PB.5275	Cleaning kit: lens fluid, lint free lens tissue/paper, brush, air blower, cotton swabs	24,00 €
PB.5276	Microscope maintenance and servicing kit, 16pcs: cleaning brush, 6 pcs screwdriver set, air blower, 3 pcs Allen key, 1.5, 2, 2.5 mm, lens cleaning fluid 20 ml, cleaning cloth 140 x 140 mm, 100 pcs Lens tissue sheets, tube of maintenance grease, 10 ml bottle of oil, packed in a nice toolbox	33,00 €
Ring LED light for NexiusZoom		
LE.1974	 Ring illuminator with 72 LEDs with adjustable light intensity. External power supply 110-240 Vac / 12 Vdc (50/60 Hz). With segment controller. Brightness of 20,000 Lux at 100 mm distance and color temperature of 6500 K. Mounting for stereo head with diameter from 25 to 61 mm	159,00 €
LE.1973	 Ring illuminator with 144 LEDs with adjustable light intensity. External power supply 100-240 Vac / 12 Vdc (50/60 Hz). With segment controller. Brightness of 23,000 Lux at 100 mm distance and color temperature of 6500 K. Mounting for stereo head with diameter from 25 to 61 mm	270,00 €
LE.1975	 Polarizing ring illuminator with 16 High Power white LEDs with adjustable light intensity. External power supply 100-240 Vac / 12 Vdc 2A 24W (50/60 Hz). Brightness of 44.000 Lux at 100 mm distance and color temperature of 6000-6500 K. Mounting for stereo head with diameter from 35 to 62 mm	335,00 €

Product	Description	EUR
DZ series		
DZ.1600	DZ stereo zoom 1:6.3 microscope with 0 to 35° ergo binocular head with EWF 10x/22 eyepieces, 1x objective, magnifications from 8x to 50x and large ergonomic LED stand  A pair of secured wide field eyepieces EWF10x/22mm with eyecups Binocular ergonomically head with 0 to 35° inclined tubes. Diopter adjustment on two tubes Interpupillary distance adjustable between 55 and 75 mm 0.8x to 5x zoom objective, magnification from 8x to 50x. Field of view from 27.5 mm to 4.4 mm. Working distance 78 mm Ergonomically stand with coarse and fine adjustments, stage plate and 2 object clamps Incident and transmitted 3W LED illumination with internal power supply Both illumination intensities can be adjusted separately Supplied with dust cover in Styrofoam case	3.779,00 €
DZ.1605	DZ stereo zoom 1:6.3 microscope with 20° fixed binocular head with EWF 10x/22 eyepieces, 1x objective, magnification from 8x to 50x and large ergonomic LED stand  A pair of secured wide field eyepieces EWF10x/22mm with eyecups Binocular head with fixed 45° inclined tubes Diopter adjustment on two tubes Interpupillary distance adjustable between 55 and 75 mm 0.8x to 5x zoom objective, magnification from 8x to 50x. Field of view from 27.5 mm to 4.4 mm. Working distance 78 mm Ergonomically stand with coarse and fine adjustments, stage plate and 2 object clamps Incident and transmitted 3W LED illumination with internal power supply Both illumination intensities can be adjusted separately Supplied with dust cover in Styrofoam case	2.696,00 €
DZ.1800	DZ stereo zoom 1:8 microscope with 0 to 35° ergo binocular head with EWF 10x/22 eyepieces, 1x objective, magnification from 8x to 64x and large ergonomic LED stand  A pair of secured wide field eyepieces EWF10x/22mm with eyecups Binocular ergonomically head with 0 to 35° inclined tubes. Diopter adjustment on two tubes. Interpupillary distance adjustable between 55 and 75 mm 0.8x to 6.4x zoom objective, magnification from 8x to 64x. Field of view from 27.5 mm to 3.4 mm. Working distance 78 mm Ergonomically stand with coarse and fine adjustments, stage plate and 2 object clamps Incident and transmitted 3W LED illumination with internal power supply Both illumination intensities can be adjusted separately Supplied with dust cover in Styrofoam case	4.013,00 €
DZ.1805	DZ stereo zoom 1:8 microscope with 20° fixed binocular head with EWF 10x/22 eyepieces, 1x objective, magnification from 8x to 64x and large ergonomic LED stand  A pair of secured wide field eyepieces EWF10x/22mm with eyecups Binocular head with fixed 45° inclined tubes Diopter adjustment on two tubes Interpupillary distance adjustable between 55 and 75 mm 0.8x to 6.4x zoom objective, magnification from 8x to 64x. Field of view from 27.5 mm to 3.4 mm. Working distance 78 mm Ergonomically stand with coarse and fine adjustments, stage plate and 2 object clamps Incident and transmitted 3W LED illumination with internal power supply Both illumination intensities can be adjusted separately Supplied with dust cover in Styrofoam case	2.930,00 €

Product	Description	EUR
DZ.1100	DZ stereo zoom 1:10 microscope with 0 to 35° ergo binocular head with EWF 10x/22 eyepieces, 1x objective, magnification from 8x to 80x and large ergonomic LED stand	5.217,00 €
	A pair of secured wide field eyepieces EWF10x/22mm with eyecups Binocular ergonomically head with 0 to 35° inclined tubes. Diopter adjustment on two tubes Interpupillary distance adjustable between 55 and 75 mm 0.8x to 8x zoom objective, magnification from 8x to 80x. Field of view from 27.5 mm to 2.8 mm. Working distance 78 mm Ergonomically stand with coarse and fine adjustments, stage plate and 2 object clamps Incident and transmitted 3W LED illumination with internal power supply Both illumination intensities can be adjusted separately Supplied with dust cover in Styrofoam case	
DZ.1105	DZ stereo zoom 1:10 microscope with 20° fixed binocular head with EWF 10x/22 eyepieces, 1x objective; magnification from 8x to 80x and large ergonomic LED stand	4.133,00 €
	A pair of secured wide field eyepieces EWF10x/22mm with eyecups Binocular head with fixed 45° inclined tubes Diopter adjustment on two tubes Interpupillary distance adjustable between 55 and 75 mm 0.8x to 8x zoom objective, magnification from 8x to 80x. Field of view from 27.5 mm to 2.8 mm. Working distance 78 mm Ergonomically stand with coarse and fine adjustments, stage plate and 2 object clamps Incident and transmitted 3W LED illumination with internal power supply Both illumination intensities can be adjusted separately Supplied with dust cover in Styrofoam case	
DZ series multi-heads		
DZ.1604	DZ stereo zoom 1:6.3 microscope with two 20° fixed binocular heads with EWF 10x/22 eyepieces, 1x objective, magnification from 8x to 50x and large ergonomic stand without illumination	10.066,00 € BTO
	A pair of secured wide field eyepieces EWF10x/22mm with eyecups Two binocular heads with fixed 45° inclined tubes Diopter adjustment on two tubes Interpupillary distance adjustable between 55 and 75 mm 0.8x to 5x zoom objective, magnification from 8x to 50x. Field of view from 27.5 mm to 4.4 mm. Working distance 78 mm Built-in LED pointer Ergonomically 350 (L) x 420 (D) mm stand with coarse and fine adjustments, 250 x 250 mm glass stage plate. Total arm length 750 mm. Without illumination Supplied with dust cover in Styrofoam case	

Product	Description	EUR	
DZ.1804	DZ stereo zoom 1:8 microscope with two 20° fixed binocular heads with EWF 10x/22 eyepieces, 1x objective, magnification from 8x to 64x and large ergonomic stand without illumination	10.687,00 €	BTO
	A pair of secured wide field eyepieces EWF10x/22mm with eyecups		
	Two binocular heads with fixed 45° inclined tubes		
	Diopter adjustment on two tubes		
	Interpupillary distance adjustable between 55 and 75 mm		
	0.8x to 64x zoom objective, magnification from 8x to 64x. Field of view from 27.5 mm to 3.4 mm. Working distance 78 mm		
	Built-in LED pointer		
	Ergonomically 350 (L) x 420 (D) mm stand with coarse and fine adjustments, 250 x 250 mm glass stage plate. Total arm length 750 mm. Without illumination		
	Supplied with dust cover in Styrofoam case		
DZ.1104	DZ stereo zoom 1:10 microscope with two 20° fixed binocular heads with EWF 10x/22 eyepieces, 1x objective; magnification from 8x to 80x and large ergonomic stand without illumination	11.343,00 €	BTO
	A pair of secured wide field eyepieces EWF10x/22mm with eyecups		
	Two binocular heads with fixed 45° inclined tubes		
	Diopter adjustment on two tubes		
	Interpupillary distance adjustable between 55 and 75 mm		
	0.8x to 8x zoom objective, magnification from 8x to 80x. Field of view from 27.5 mm to 2.8 mm. Working distance 78 mm		
	Built-in LED pointer		
	Ergonomically 350 (L) x 420 (D) mm stand with coarse and fine adjustments, 250 x 250 mm glass stage plate		
	Total arm length 750 mm		
	Without illumination		
	Supplied with dust cover in Styrofoam case		
DZ series - Individual components			
DZ.0630	DZ Zoom 1:6.3 body without head, eyepieces, objectives and stand	1.187,00 €	
DZ.0800	DZ Zoom 1:8 body without head, eyepieces, objectives and stand	1.420,00 €	
DZ.1000	DZ Zoom 1:10 body without head, eyepieces, objectives and stand	2.624,00 €	
DZ.2020	DZ Binocular ergonomically head with inclined tubes from 0-35° without eyepieces	1.487,00 €	
DZ.2025	DZ Binocular head with 20° inclined tubes without eyepieces	405,00 €	
DZ.3010	Eyepiece EWF 10x/22 mm (1pc) for DZ-series head	118,00 €	
DZ.3012	Eyepiece EWF 10x/22 mm with crosshair with 10/100 micrometer scale for DZ-series head (1pc)	161,00 €	
DZ.3015	Eyepiece EWF 15x/16 mm (1pc) for DZ-series head	205,00 €	
DZ.3020	Eyepiece EWF 20x/12 mm (1pc) for DZ-series head	212,00 €	
DZ.3025	Pair of eyecups for DZ series	13,30 €	
DZ.4003	Plan achromatic common main 0.3x objective for DZ-series. Working distance 276 mm	378,00 €	
DZ.4005	Plan apochromatic 0.5x common main objective for DZ-series. Working distance 126 mm	771,00 €	
DZ.4010	Plan apochromatic 1x common main objective for DZ-series. Working distance 78 mm	285,00 €	
DZ.4020	Plan apochromatic 2x common main objective for DZ-series. Working distance 33 mm	2.004,00 €	
DZ.5020	Universal stand with holder without illumination for DZ-series	513,00 €	
DZ.5040	Large Ergonomic LED stand with head holder for DZ-series	690,00 €	
DZ.9005	Single photo tube attachment for DZ-series with one 23.2 mm tube	516,00 €	
DZ.9007	Dual photo tube attachment for DZ-series with two 23.2 mm tube	1.119,00 €	
AE.5128	SLR camera adapter with 2x lens for 23.2 mm tubes of DZ-series	335,00 €	
DZ.9010	Dual Iris diaphragm for DZ-series	259,00 €	
DZ.9011	C-mount adapter with 0.4x lens	287,00 €	

Product	Description	EUR
DZ.9013	C-mount adapter with 1x lens	75,00 €
DZ.9020	X-Y stage 185 x 145 mm for DZ-series with 50 x 50 mm travel	255,00 €
DZ.9030	Polarization set for DZ-series	193,00 €
DZ.9040	Darkfield attachment for DZ.5040	69,00 €
DZ.9045	Analyzer 360° rotatable, fits the DZ-series objectives	98,00 €
DZ.9047	Polarization filter 360° rotatable, fits on DZ.5040 stand	131,00 €
DZ series - Fluorescence attachment		
DZ.9050	Fluorescence attachment with GFP-B and GFP-L filter sets and one 23.2 mm photo tube	10.260,00 €
DZ.9060	Filter set GFP-B (EX460-500, DM505, EM510-560)	605,00 €
DZ.9061	Filter set GFP-L (EX460-500, DM505, EM510LP)	873,00 €
DZ.9749	Empty filter block for DZ series fluorescence attachment	199,00 €
DZ series - Accessories		
AE.1900	Glass object-plate, ø 100mm	38,00 €
AE.1903	Black/White stage plate, ø 100mm	38,00 €
SL.1385	100 W Mercury lamp	312,00 €
AE.1385	100 W Mercury lamp	312,00 €
SL.5520	LED replacement for incident illumination for DZ.5040	41,00 €
AE.5130	Universal SLR adapter with built-in 2x lens for standard 23.2 mm tube. Needs T2 adapter	261,00 €
AE.5025	T2 ring for Nikon D SLR digital camera	38,00 €
AE.5040	T2 ring for Canon EOS SLR digital camera	38,00 €
AE.1112	Micrometer 76 x 26 mm slide, 50 mm/500 divisions, 100 µm per division	76,00 €
LE.1974	Ring illuminator with 72 LEDs with adjustable light intensity. External power supply 110-240 Vac / 12 Vdc (50/60 Hz). With segment controller. Brightness of 21,000 Lux at 100 mm distance and color temperature of 6500 K. Mounting for stereo head with diameter from 25 to 61 mm	159,00 €
		
LE.1973	Ring illuminator with 144 LEDs with adjustable light intensity. External power supply 100-240 Vac / 12 Vdc (50/60 Hz). With segment controller. Brightness of 23,000 Lux at 100 mm distance and color temperature of 6500 K. Mounting for stereo head with diameter from 25 to 61 mm	270,00 €
		
LE.1975	Polarizing ring illuminator with 16 High Power white LEDs with adjustable light intensity. External power supply 100-240 Vac / 24 Vdc 220 mA 4.5W (50/60 Hz). Brightness of 44.000 Lux at 100 mm distance and color temperature of 6000-6500 K. Mounting for stereo head with diameter from 35 to 62 mm	335,00 €
		
PB.5245	Lens cleaning paper, 100 sheets per pack	3,20 €
PB.5274	Iso propyl alcohol 99%, 200 ml	29,00 €
PB.5275	Cleaning kit: lens fluid, lint free lens tissue/paper, brush, air blower, cotton swabs	24,00 €